Microsart® maxi.vac | Microsart® mini.vac

Laboratory Vacuum Pumps for Microbiology Analysis

Vacuum pumps

Microsart® maxi.vac Laboratory Vacuum pump 22 l/min for 3- and 6-place manifolds

Description

Microsart® maxi.vac and Microsart® mini.vac are especially designed to meet the requirements of microbiology laboratories. The absolute vacuum is limited to 100 mbar. The adjustable vacuum allows for trouble free working according to ISO 8199 ("Water quality: General guide to enumeration of micro-organisms by culture). The differential pressure during filtration should not exceed 700 mbar to allow a gentle enrichment of the microorganisms.

Microsart® mini.vac Laboratory Vacuum pump 6 l/min for 1-branch manifolds and individual filter holders

Efficient filtration

The intelligent pump design facilitates easy and efficient operation, even in limited working space. The noise development during operation is reduced by 50% compared to conventional membrane pumps. An integral part of the pumps is a newly structured membrane which guarantees leak-proof operation.

The pumps are available in two sizes. The small compact pump Microsart® mini.vac with 6 l/min air flow is ideal for single vacuum filtrations and the standard pump Microsart® maxi.vac, with 22 l/min air flow, for 3- to 6-place manifolds.

- Compliance with ISO 8199 "Water Testing Guideline"
- Adjustable vacuum
- Reduced noise
- Structured Membrane – Patented protection against leakage
- Compact design – Small foot print
- Maintenance-free, 5,000 operating hours guaranteed

Technical Specifications


Microsart® maxi.vac Laboratory Vacuum Pump 22 l

Delivery	22 (l/min)
Absolute Vacuum	100 (mbar)
Noise level [100 mbar]	57.5 – 59.0 dBA
Operating Pressure (max)	1 (bar)
Materials (contact with filtrate possible)	Aluminum, CR (Neoprene), NBR (Perbunan)
Connectors for Tube	ID 9 (mm)
Ambient Temperature	5 ... 40°C
Mains 16694-2-50...	230 V/50 Hz
Mains 16694-1-60...	115 V/60 Hz
Motor Protection	IP 44
Power P1	130 W
Operating Current	0.9 A
Weight	7.1 kg
Dimensions L×H×W (mm)	261×204×110

Microsart® mini.vac Laboratory Vacuum Pump 6 l

Delivery	6 (l/min)
Absolute Vacuum	100 (mbar)
Noise level [100 mbar]	53.5 dBA
Operating Pressure (max)	2.5 (bar)
Materials (contact with filtrate possible)	PPS, EPDM, FPM (Viton)
Connectors for Tube	ID 4 (mm)
Ambient Temperature	5 ... 40°C
Mains 16694-2-50...	230 V/50 Hz
Mains 16694-1-60...	115 V/60 Hz
Motor Protection	IP 20
Power P1	65 W
Operating Current	0.63 A
Weight	1.9 kg
Dimensions L×H×W (mm)	164×141×90

Pump capacity at stable vacuum 100 mbar

The graph illustrates the maximum potential of the Laboratory Vacuum pumps during filtration with standard Cellulose Nitrate membranes of different pore sizes. In example, the 6 l vacuum pump delivers enough air flow for three filtrations in parallel with 0.45 µm membranes.

Ordering Information Microsart® maxi.vac Laboratory Vacuum Pump 22 l

Order Number	Mains	Delivery
16694-2-50-22	230 V	22 l/min
16694-1-60-22	115 V	22 l/min


Spare Parts Microsart® maxi.vac

Order Number	Description
1EV---0002	Fine adjustment head incl. vacuometer
1EH---0002	Sound absorber
1ED---0055	Repair part set incl. membrane, sealing ring, valves

Ordering Information Microsart® mini.vac Laboratory Vacuum Pump 6 l

Order Number	Mains	Delivery
16694-2-50-06	230 V	6 l/min
16694-1-60-06	115 V	6 l/min

Characteristic set-up of a laboratory vacuum filtration

More information about vacuum filtration equipment is available at:
www.sartorius-stedim.com/microbio – Colony Counting – Combisart Hardware Set-Ups

Spare Parts Microsart® mini.vac

Order Number	Description
1EV---0001	Fine adjustment head incl. vacuometer
1EH---0001	Sound absorber
1ED---0054	Repair part set incl. membrane, sealing ring, valves
1EV---0003	Fine adjustment head incl. manometer

Microsart® e.jet

Laboratory Transfer Pump with Quick Connection for Microbiological Analysis

Typical Set-up: Microsart® e.jet connected to a 3-branch manifold with Microsart® Funnels 100

New pump for vacuum filtration in Microbiological Analysis

The Microsart® e.jet is a laboratory vacuum pump which creates vacuum and concomitantly transfers the filtered liquid to waste. The Microsart® e.jet is ideal for sample preparation in microbiological analysis achieving a trans membrane pressure of 600 mbar and a flow rate of 4.0 NI/min (4.0 Normliters water displacement by air in one minute). Constant flow rates and a defined maximum vacuum guarantees a smooth and reliable filtration.

Reduced operating complexity

Until now vacuum equipment for the Membrane Filtration Method consists of numerous parts including connectors, tubes, vacuum containers, protection filter, Woulff's bottle and a vacuum pump. After several samples the vacuum must be broken to empty the filtrate collection container. The complete traditional equipment requires far more laboratory space and is time consuming to operate and maintain. The new Microsart® e.jet greatly reduces operating complexity.

Smart design fits any laboratory environment

The Microsart® e.jet pump is an ideal accessory for 3-branch and 1-branch manifolds. Compared to traditional equipment, the

Microsart® e. jet and stainless steel manifold require only 30% of the average space, which translates to less congestion when working in laminar flow cabinets.

Safe operation | Maintenance-free

Traditional vacuum pumps often lose their efficiency and capability to generate sufficient vacuum, when liquid is drawn into the pump head. The Microsart® e.jet is designed to pump both gas and liquids, meaning no loss of efficiency or malfunction from water drawn into the pump head.

Building-up the vacuum filtration system is easy and fast thanks to the innovative Quick Connections. The Microsart® e.jet Transfer Pump is equipped with Quick Connection Nipples assembled to Quick Connection Couplings on hose nipples for DN 10 tubings. Simply push-to-connect for assembling and pull-to-disassembling the whole system within seconds. The Quick Connections are non-shut-off.

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