

Portable Incubator

◎ Can carry the training equipment everywhere.



Patent No. : ZL 2018 3 0095051.6

Main Technical Parameters:

Model		DH2500AB
Cycle Mode		Natural convection
Function	Tem. Range	RT+5-50°C
	Tem. Resolution Ratio	0.1°C
	Tem. Motion	± 0.5°C
	Tem. Uniformity	± 1.5°C
Structure	Inner Chamber	Mirror stainless steel
	Outer Shell	Cold rolling steel electrostatic spraying exterior
	Insulation layer	High quality foam board
	Heater	Heating Wire
	Power Rating	0.3KW
Controller	Tem. control mode	PID Intelligent
	Tem. setting mode	Touch button setting
	Tem. display mode	Measuring temperature: LED upper row ; Setting temperature: the lower row
	Timer	0-9999min (with timing wait function)
	Operation function	Fixed temperature operation, timing function, auto stop.
	Additional function	Sensor deviation correction, temperature overshoot self-tuning, internal parameter locking, power-off parameter memory
	Sensor	CU50
Safety device		Over temperature sound-light alarm
Specification	Inner Chamber size(W*L*H)(mm)	200*160*230
	Exterior size (W*L*H)(mm)	300*330*330
	Packing size (W*L*H)(mm)	340*370*390
	Volume	7.4L
	Supply(50/60HZ)Current rating	AC220V/1.3A
	NW/GW (kg)	9/11
Optional accessories		Vehicle supply/DC12V supply/more smaller size choice