

Secador® Desiccator Cabinets



Durable, co-polyester construction blocks 99% of UV light (mini is polystyrene)

Patented door seal provides superior protection from moisture

Portable & stackable



NOTE: All 1.0, 2.0, 3.0 and 4.0 cabinets are interchangeably stackable, except 4.0 horizontal.

Built-in, convenient carrying handles

Easy-to-view dial hygrometer in door

Automated models provide continuous desiccant regeneration*



Door can be secured with tamper-evident tags or padlock

Automated models available in 100VAC, 120VAC and 230VAC

All sizes are available with optional gas ports to allow gas purging of the cabinet.

* All non-electric models can be used with any desiccant including Disposable Desiccant Cartridges F42048-0065, -0100, or Desi Can™ Reusable Desiccant Canister F42046-0000.

Secador®

1.0, 2.0, 3.0 and 4.0 Vertical Desiccator Cabinets

Interchangeably Stackable up to 3 Units High

U.S. Patents: 6,834,920; 7,318,630 and D474,626



- All units have a compact 34.1W x 41.4cmD (13.4 x 16.3") footprint
- Doors open to the left
- Different colors are available to suit your needs

Size Specifications

Model	Internal Volume	Holds (12oz. soda cans)	# of Shelves	Height
1.0	0.74 cu. ft.	20	1 (fixed)	21.3cm (8.4")
2.0	1.16 cu. ft.	40	2	31.4cm (12.4")
3.0	1.55 cu. ft.	60	3	41.7cm (16.4")
4.0	1.87 cu. ft.	70	3	51.8cm (20.4")

Secador® 1.0, 2.0, 3.0 and 4.0 Desiccator Cabinets

CATALOG NO.	MODEL	DESCRIPTION	PRICE/EACH
F42071-0000	1.0	Clear	\$208.40
F42071-0006	1.0	Clear w/Blue End-caps	208.40
F42071-0007	1.0	Blue w/Clear Door	208.40
F42071-0008	1.0	Amber (50+% visible light blocking.)	234.70
F42072-1000	2.0	Clear	416.20
F42072-1008	2.0	Amber (50+% visible light blocking.)	437.40
F42073-1000	3.0	Clear	530.10
F42073-1008	3.0	Amber (50+% visible light blocking.)	548.70
F42074-1000	4.0	Clear	737.55
F42074-1006	4.0	Clear, with Blue End-caps	737.55
F42074-1007	4.0	Blue with Clear Door	737.55
F42074-1008	4.0	Amber (50+% visible light blocking.)	751.30

Secador® Cabinet Replacement Shelves p. 109

For more detailed information on the full Secador® line see page 88.